



Introducing the new Testogel 16.2 mg/g gel

The number one selling testosterone replacement therapy worldwide¹ now comes in an easy-to-use pump.²

Compared to Testogel 50 mg³, gel in sachet (10mg/g), Testogel 16.2 mg/g gel:

- Offers the ability to titrate more easily and accurately
- Delivers a smaller, more concentrated volume of gel, enabling patients to apply more easily⁴

Find out more by visiting
www.testogelpump.co.uk



*Put back the testosterone,
complete the man*

Adverse events should be reported. Reporting forms and information can be found at www.mhra.gov.uk/yellowcard Adverse events should also be reported to Besins Healthcare (UK) Ltd, 28 Poland Street, London, W1F 8QN. Tel: 0203 862 0920. Email: pharmacovigilance@besins-healthcare.com

References: 1. IMS Health. Internal calculations based on IMS Health, IMS MIDAS MAT Q3 2017 (LCD/MNF). 2. Testogel 16.2 mg/g gel Summary of Product Characteristics. 3. Testogel 50 mg gel SPC 4. Kaufman JM et al. J Sex Med. 2012; 9(4): 1149-1161. Date of Preparation: January 2018. TES/2018/004.

Full prescribing information is available at this stand